

ACCEPTED
FOR REGISTRATION

December 23, 2015

DOC ID 543986

New York State Department
of Environmental Conservation
Division of Materials Management
Pesticide Product Registration

microdot[®]

MINUTE WIPE

ONE STEP CLEAN and DISINFECT 1 MINUTE

Virucidal* · Bactericidal · Tuberculocidal · Fungicidal

1 MINUTE STEP

ACTIVE INGREDIENTS:
Ethyl Alcohol 72.50%
Didecyl Dimethyl Ammonium Chloride .33%

OTHER INGREDIENTS*
Total 27.17%
100.00%

* Other ingredients do not include the weight of the dry wipes.

KEEP OUT OF REACH OF CHILDREN
CAUTION

NET CONTENTS:
1 lb. 8 ounces
160 Pre-saturated 6 x 6.75 in
WIPES

EPA Reg. No. 88494-2-88459 (See side panel for additional precautionary statements)
EPA Est. No. 85743-TUR-1

ONE MINUTE • ONE STEP

Disinfecting Wipe for Hard, Non-Porous Surfaces

Patient Care Equipment and Point-of-Care Equipment

Kills Viruses* Most Associated with the Common Cold • Kills Antibiotic Resistant Bacteria MRSA and VRE

Kills Flu Virus* • Kills Pandemic 2009 H1N1 Influenza A virus*

UNSCENTED FRAGRANCE FREE BLEACH-FREE

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

DISPENSING INSTRUCTIONS: To Dispense Wipes: Open lid, remove entire cover, remove and discard seal. From center of wipes roll pull up a corner of the wipe, twist it and thread it through the hole in the container cover. Pull through about one inch. Replace cover tightly. Pull out first wipe and tear off at a 45° angle. Remaining wipes feed automatically, ready for next use. When not in use, keep the tip of the exposed wipe inside the well in the cover, and the lid snapped closed.

CLEANING: Use one or more wipes, as necessary, to wet surface sufficiently and to thoroughly clean the surface.

DISINFECTING: To disinfect hard, non-porous surfaces, use one or more wipes, as necessary, to thoroughly wet the surface to be treated. Treated surface must remain visibly wet for 1 minute to achieve complete disinfection of all pathogens listed on this label. For visibly soiled surfaces, clean before following disinfecting instructions. Allow surfaces to air dry. For surfaces that may come in contact with food, a potable water rinse is required.

***VIRUCIDAL:** microdot® Minute Wipe is an effective virucide on hard, non-porous surfaces, in the presence of an organic load, against Poliovirus serotype 1 [PV-1], Hepatitis B Virus, Hepatitis C Virus, Rhinovirus [ATCC VR-340], Norovirus, Human Coronavirus [ATCC VR-1558], Rotavirus Strain WA, Respiratory Syncytial Virus [ATCC VR-26], Human Immunodeficiency Virus [HIV-1], Herpes Simplex Virus Type 1 [ATCC VR-539], and Influenza Virus type A2 [ATCC VR-544], when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

KILLS HIV-1, HBV and HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 [HIV-1], Hepatitis B Virus [HBV], and Hepatitis C Virus [HCV]

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV, and HCV FOR SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS.

PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks and eye coverings. **CLEANING PROCEDURES:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of microdot® Minute Wipe. **CONTACT TIME:** Allow surface to remain wet for 1 full minute. **DISPOSAL OF INFECTIOUS MATERIALS:** Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

BACTERICIDAL: microdot® Minute Wipes an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against E. coli, Pseudomonas aeruginosa, Salmonella enterica, Staphylococcus aureus, Streptococcus pyogenes, Vancomycin resistant Enterococcus faecium, Methicillin resistant Staphylococcus aureus, Klebsiella pneumoniae, Listeria monocytogenes, Acinetobacter baumannii, Campylobacter jejuni, and Burkholderia cepacia, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry. **TUBERCULOCIDAL:** microdot® Minute Wipe is an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against Mycobacterium bovis, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry. **FUNGICIDAL:** microdot® Minute Wipes an effective fungicide against Trichophyton mentagrophytes, when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, exercise facilities or other clean, hard, non-porous, surfaces commonly contacted by bare feet. The surface must remain visibly wet for 1 minute. Allow surfaces to air dry. **EFFICACY TESTS HAVE DEMONSTRATED THAT microdot® Minute Wipe IS AN EFFECTIVE, VIRUCIDE, BACTERICIDE, TUBERCULOCIDE AND FUNGICIDE IN THE PRESENCE OF ORGANIC SOIL**

FOOD PROCESSING AND FOOD PACKAGING FACILITIES: Before use in federally inspected meat and poultry food processing plants and dairies, food products and packaging materials must be removed from the room or carefully protected. microdot® Minute Wipe will not harm most surfaces: including stainless steel, aluminum, plastics, fiberglass, vinyl, and most painted surfaces. For plastic and painted surfaces, spot test on an inconspicuous area before use. Do not use on unpainted wood. **Not for use on people or pets.**

PRECAUTIONARY STATEMENTS - Hazards to Humans and Domestic Animals:

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet.

FIRST AID

IF IN EYES: Hold eyes open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center, doctor, or going for treatment. For emergency information concerning this product, call the **National Pesticides Information Center at 1-800-858-7378** seven days a week, 6:30 am to 4:30 pm Pacific Time.

PHYSICAL OR CHEMICAL HAZARDS: Flammable. Keep away from heat or open flame.

STORAGE AND DISPOSAL - Do not contaminate water, food or feed by storage or disposal.

PESTICIDE STORAGE: Always store this product and all pesticides in the original container. If a leaky container must be contained within another, mark the outer container to identify the contents. Protect from extreme heat or cold. **PESTICIDE DISPOSAL:** Waste resulting from the use of this product must be disposed of on site or at an approved waste disposal facility. **CONTAINER DISPOSAL:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available, or discard empty container in trash..

CS
CAMBRIDGE SENSORS USA, LLC
©2015 Cambridge Sensors USA, LLC.

Manufactured in Turkey for
Cambridge Sensors USA, LLC.
10051 Bode Road Plainfield, Illinois 60585
info@microdotcs.com • www.microdotcs.com
Customer Service: Toll Free 877.374.4062 (M-F 9am-5pm CST)

