

February 26, 2018

DOC ID 555769

New York State Department
of Environmental Conservation
Division of Materials Management

Pesticide Product Registration

PERFORMANCE WIPES™

KILLS 99.9% OF BACTERIA IN JUST 15 SECONDS**!

Antibacterial Formula Cleans and Disinfects

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling.

FIRST AID

Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

IF IN EYES: Hold eyes open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

CONTAINER STORAGE/DISPOSAL. Store in an area inaccessible to children. **Wipes:** Dispose of wipe in trash after use. Do not flush. **Dispenser:** Refill only with Performance Wipes. Do not reuse or refill except as described in the directions for use. Keep dispenser closed to retain moisture. **Plastic Bag:** Do not reuse or refill empty plastic bag. Offer empty plastic bag for recycling if available or discard in trash.

** Staphylococcus aureus (Staph), Klebsiella pneumoniae, Escherichia coli (E.coli) O157:H7

Manufactured for:



19 Bennett Street
Lynn, MA 01905
800-225-9473
781-593-4000

www.ercwipe.com

EPA REG. No. 1839-190-85282
EPA EST. No. 87871-NY-001

Net Contents: 11 lbs., 0 oz.
(800 Pre-moistened Wipes (6" x 6"))

- Kills Cold¹ and Flu Viruses†††
- Kills HIV-1
- Kills Antibiotic Resistant Bacteria:

CA-MRSA, MRSA, VISA, VRE and NDM-1

- Alcohol, Solvent, and Bleach Free

††† Influenza A Virus, Pandemic 2009 H1N1
Influenza A Virus, Human Coronavirus¹

Active Ingredients

Alkyl (60% C ₁₄ , 30% C ₁₆ , 5% C ₁₂ , 5% C ₁₈)	dimethyl benzy ammonium chlorides.....	0.14%
Alkyl (68% C ₁₂ , 32% C ₁₄) dimethyl ethylbenzy ammonium chlorides.....		0.14%
Inert Ingredients		99.72%
Total		100.00%

**KEEP OUT OF REACH
OF CHILDREN**

CAUTION

See side panels for additional precautionary statements.

Performance Wipes are non-woven, disposable, wipes, pre-saturated with a quaternary disinfectant cleaner. This product is intended for use in hospitals and healthcare settings, critical care areas, veterinary clinics, laboratory, transport vehicles, health care facilities, shops, gyms, schools, shopping centers, catering facilities, and daycare centers. This product is formulated to disinfect hard, non-porous, inanimate environmental surfaces such as floors, walls, metal surfaces, stainless steel surfaces, glazed porcelain, glazed ceramic tile, plastic surfaces, bathrooms, shower stalls, bathtubs, and cabinets.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. This product is not for use on critical and semi-critical medical device surfaces.

CLEANING: Use one Performance Wipe to completely remove all loose debris from surfaces to be treated.

DISINFECTING: To disinfect hard, non-porous surfaces, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for disinfecting. Use enough wipes for treated surface to remain visibly wet for 5 minutes. For heavily soiled surfaces, clean before following disinfecting instructions. For surfaces to air dry. Performance Wipes are an effective disinfectant against Escherichia coli O157:H7 (ATCC 4388), Staphylococcus aureus (ATCC 13482), Salmonella enterica (ATCC 10708), Staphylococcus aureus (ATCC 6538), Streptococcus pyogenes (Strp) (ATCC 19619), Vancomycin resistant Enterococcus faecium (VRE) (ATCC 52059), Methicillin resistant Staphylococcus aureus (MRSA) (ATCC 33591), Community Associated Methicillin resistant Staphylococcus aureus (CA-MRSA), Genotype USA400, Community Associated Methicillin resistant Staphylococcus aureus (CA-MRSA) genotype USA300 and Vancomycin intermediate resistant Staphylococcus aureus (USA) (ATCC 70629) and *Klebsiella pneumoniae* New Delhi, Metallo-Beta Lactamase (NDM-1) Capampium Resistant (CDC Clinical Isolate No. 1002).

***VIRUCIDAL ACTIVITY:** Performance Wipes are an effective virucide on hard, non-porous surfaces against Influenza A2/Japan, Herpes Simplex Virus 1 (ATCC VR-260), Human Immunodeficiency Virus, Type 1 (HIV-1) and Herpes Simplex Virus 2 (ATCC VR-234), Duck Hepatitis B Virus (DHBV), Human Hepatitis B Virus (HBV), Bovine Viral Diarrhea Virus (BVD) and Human Hepatitis C Virus (HCV) when the treated surface is allowed to remain wet for 5 minutes. Allow surface to air dry. Performance Wipes is an effective virucide on hard, non-porous surfaces against Human Coronavirus (ATCC VR-740) when the treated surface is allowed to remain wet for 10 minutes. Allow surface to air dry.

KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings (Hospital, Nursing Home) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects may be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1, Human Hepatitis B Virus (HBV) and Human Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV AND HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks, and eye coverings.

CLEANING PROCEDURES: Blood and other body fluids must be thoroughly cleaned from surfaces, and objects before application of this product.

CONTACT TIME: Allow surface to remain wet for 5 minutes.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

FUNGICIDAL ACTIVITY: Performance Wipes are an effective fungicide against Trichophyton mentagrophytes (ATCC 9633) (causative agent of Athlete's Foot Fungus), when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, exercise facilities or other clean, hard, non-porous, surfaces commonly contacted by bare feet. The surface must remain visibly wet for 3 minutes. Allow surface to air dry.

EFFICACY TESTS HAVE DEMONSTRATED THAT PERFORMANCE WIPES ARE AN EFFECTIVE BACTERICIDE, VIRUCIDE, AND FUNGICIDE IN THE PRESENCE OF ORGANIC SOIL (6% BLOOD SERUM).

SANITIZING/DEODORIZING: To sanitize/deodorize hard, non-porous surfaces, wipe surface, using enough wipes to allow treated surface to remain wet for 15 seconds for Staphylococcus aureus (Staph) (ATCC 6538), Klebsiella pneumoniae (ATCC 4382), Listeria monocytogenes (ATCC 35182), Escherichia coli (E. coli) O157:H7 (ATCC 43885) and 1 minute for Campylobacter jejuni (ATCC 29429). Let surface dry. For highly soiled surfaces, clean excess dirt first. For surfaces that may come in contact with food, a potable water rinse is required. For plastic and painted surfaces, spot test on an inconspicuous area before use. Performance Wipes will not treat most surfaces: including acrylic, fiberglass and vinyl. Do not use on unpainted wood.

Not For Personal Cleansing.

This is NOT A Baby Wipe!

Do Not Flush in Toilet.

TO OPEN PACKAGE: Cut a "Y" out of the top of the bag. Pull out the first wipe approximately 6 to 8 inches from the center of the roll. Place entire bag inside of center pull dispenser. Thread the wipe through the opening on the dispenser. Close the lid of the dispenser. Dispenser may be mounted with wipes dispensing up.

