

# Super

# Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

EXTRA LARGE  
WIPE

ACCEPTED  
FOR REGISTRATION

June 8, 2015

New York State Department  
of Environmental Conservation  
Division of Materials Management  
Pesticide Product Registration

Doc id: 542414



## BACTERICIDAL, TUBERCULOCIDAL, AND VIRUCIDAL\* IN 2 MINUTES

### ACTIVE INGREDIENTS

n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chlorides .....	0.25%
n-Alkyl (60% C <sub>12</sub> , 30% C <sub>14</sub> , 5% C <sub>12</sub> , 5% C <sub>14</sub> ) dimethyl benzyl ammonium chlorides ....	0.25%
Isopropyl Alcohol .....	55.00%

OTHER INGREDIENTS .....	44.50%
TOTAL .....	100.00%

Does not include the weight of the wipe.

KEEP OUT OF REACH OF CHILDREN

**WARNING**

See back panel for additional  
precautionary statements

**CONTAINS: 65 Wipes** 7.5 x 15 in (19.0 x 38.1 cm)  
Net Wt. 2 lb 1 oz (938g)

**+PDI**

REORDER NO. Q86984

## Effective against:

**BACTERIA** • *Bordetella pertussis* • *Burkholderia cepacia*  
• *Campylobacter jejuni* • *Escherichia coli* (E.coli)  
• *Escherichia coli* O157:H7 • *Klebsiella pneumoniae*  
• *Pseudomonas aeruginosa* • *Salmonella enterica*  
• *Staphylococcus aureus*

**MULTI-DRUG RESISTANT BACTERIA** • *Acinetobacter baumannii* • *Enterobacter cloacae* NDM-1 positive  
• ESBL Producing *Escherichia coli* (E. coli) • *Klebsiella pneumoniae* - KPC-2 positive, ST258 • Methicillin Resistant *Staphylococcus aureus* (MRSA) • Vancomycin Resistant *Enterococcus faecalis* (VRE)

**VIRUSES\*** • \*Adenovirus Type 5 • \*Herpes Simplex Virus Type 2 • \*Human Coronavirus Strain 229E • \*Influenza A/Hong Kong • \*Influenza A (H1N1) virus • \*Respiratory Syncytial Virus (RSV) • \*Rhinovirus • \*Rotavirus Strain WA • \*Vaccinia virus • Kills Pandemic 2009 H1N1 Influenza A virus

**BLOODBORNE PATHOGENS** • \*Hepatitis B virus (HBV) - Duck HBV† • \*Hepatitis C virus (Human) (HCV) - Bovine Diarrhea Virus† • \*HIV-1 (AIDS virus)†

**TB** • *Mycobacterium bovis* - BCG (Tuberculosis)

**PATHOGENIC FUNGI** • *Candida albicans*

8HQ86984



(01)10310819000588

\*KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature [68°-77°F] in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

**Areas of Use: Hospital and Healthcare Settings:** Ambulatory Surgical Centers (ASC), Clinics, Dental Offices, Dialysis Clinics, Home Health Care, Hospices, Hospitals, Laboratory, Nursing homes, Physical therapy, Physician's offices, Radiology, Rehabilitation, Transport vehicles. **Critical Care Areas:** CCU, Emergency Rooms, ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU), Surgery and Surgical Intensive Care Units (SICU). **Hospital, Healthcare, and Critical Care Use Sites:** May be used on hard non-porous surfaces of: bed railings; blood glucose meters; cabinets; carts; chairs; counters; dental unit instrument trays; exam tables; gurneys; isolettes; IV poles; stethoscopes; stretchers; tables; telephones; toilet seats; and hard non-porous outside surfaces of: amalgamators and dental curing lights; diagnostic equipment; patient monitoring equipment; patient support and delivery equipment.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**DIRECTIONS FOR USE:** It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**To Dispense Wipes:** Remove lid and discard inner seal from canister. Find center of wipe roll, remove first wipe for use, twist corner of next wipe into a point and thread through the hole in the canister lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out at an angle. When not in use keep center cap of lid closed to prevent moisture loss.

**TO DISINFECT AND DEODORIZE:** To disinfect nonfood contact surfaces only: Unfold a clean wipe and thoroughly wet surface. Allow treated surface to remain wet for a full two (2) minutes. Let air dry. For heavily soiled areas, use a wipe to pre-clean prior to disinfecting. These directions also apply to *Mycobacterium bovis* BCG (Tuberculosis) at 68°F [20°C].

**\*SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS**

**Personal protection:** When using this product, wear disposable protective gloves, protective gowns, masks, and eye coverings when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids. **Cleaning procedure:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold and use first germicidal wipe to remove heavy soil. **Disposal of infectious materials:** Used wipe, blood and other body fluids should be disposed of according to local regulation for infectious waste disposal. **Contact time:** Use second germicidal wipe to thoroughly wet surface. Allow to remain wet two (2) minutes, let air dry.

**PRECAUTIONARY STATEMENTS Hazards to Humans and Domestic Animals**

**WARNING.** Causes substantial but temporary eye damage. Do not get in eyes or on clothing. Avoid contact with skin. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using restroom. Remove and wash contaminated clothing before reuse.

**FIRST AID** Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes.

**PHYSICAL OR CHEMICAL HAZARD** Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

Test wipe on small inconspicuous area first.

**STORAGE AND DISPOSAL** Do not contaminate water, food, or feed by storage and disposal. Storage: Do not store near heat or open flame. When not in use keep center cap of lid closed to prevent moisture loss. **Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Container Disposal:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling. If recycling is not available, put in trash collection.

**EPA REG. NO. 9480-4**

**EPA EST. NO. A= 9480-NY-1, C= 72956-AR-1, D= 8251-WI-4**

**Alpha character will precede batch code on product**

Manufactured for:

Professional Disposables International, Inc.

Two Nice-Pak Park

Orangeburg, NY 10962-1376

1-800-999-6423

Made in USA with domestic and imported materials

**NOT FOR USE ON SKIN FOR USE ON HARD, NON POROUS SURFACES ONLY**

# Super

# Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

ACCEPTED  
FOR REGISTRATION

June 8, 2015

New York State Department  
of Environmental Conservation  
Division of Materials Management  
Outside Product Registration

LARGE WIPE



**BACTERICIDAL, TUBERCULOCIDAL, AND VIRUCIDAL\* IN 2 MINUTES**

#### ACTIVE INGREDIENTS

n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chlorides .....	0.25%
n-Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> ) dimethyl benzyl ammonium chlorides ....	0.25%
Isopropyl Alcohol .....	55.00%
<b>OTHER INGREDIENTS</b> .....	44.50%
<b>TOTAL</b> .....	100.00%

Does not include the weight of the wipe.

**KEEP OUT OF REACH OF CHILDREN**

**WARNING** See back panel for additional  
precautionary statements

**CONTAINS: 160 Wipes** 6 x 6.75 in (15.2 x 17.1 cm)

Net Wt. 1 lb 13 oz (833g)



REORDER NO. Q55172

#### Effective against:

**BACTERIA** • *Bordetella pertussis* • *Burkholderia cepacia*  
• *Campylobacter jejuni* • *Escherichia coli* (E.coli)  
• *Escherichia coli* 0157:H7 • *Klebsiella pneumoniae*  
• *Pseudomonas aeruginosa* • *Salmonella enterica*  
• *Staphylococcus aureus*

**MULTI-DRUG RESISTANT BACTERIA** • *Acinetobacter baumannii* • *Enterobacter cloacae* NDM-1 positive  
• ESBL Producing *Escherichia coli* (E. coli) • *Klebsiella pneumoniae* - KPC-2 positive, ST258 • Methicillin Resistant *Staphylococcus aureus* (MRSA) • Vancomycin Resistant *Enterococcus faecalis* (VRE)

**VIRUSES\*** • Adenovirus Type 5 • Herpes Simplex Virus Type 2 • Human Coronavirus Strain 229E • Influenza A/Hong Kong • Influenza A (H1N1) virus • Respiratory Syncytial Virus (RSV) • Rhinovirus • Rotavirus Strain WA • Vaccinia virus • Kills Pandemic 2009 H1N1 Influenza A virus

**BLOODBORNE PATHOGENS** • Hepatitis B virus (HBV) - Duck HBV† • Hepatitis C virus (Human) (HCV) - Bovine Diarrhea Virus†  
• HIV-1 (AIDS virus)†

**TB** • *Mycobacterium bovis* - BCG (Tuberculosis)

**PATHOGENIC FUNGI** • *Candida albicans*

8LQ55100

\*KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature (68°-77°F) in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids; and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

**Areas of Use: Hospital and Healthcare Settings:** Ambulatory Surgical Centers (ASC), Clinics, Dental Offices, Dialysis Clinics, Home Health Care, Hospices, Hospitals, Laboratory, Nursing homes, Physical Therapy, Physician's offices, Radiology, Rehabilitation, Transport vehicles. **Critical Care Areas:** CCU, Emergency Rooms, ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU), Surgery and Surgical Intensive Care Unit (SICU). **Hospital, Healthcare, and Critical Care Use Sites:** May be used on hard non-porous surfaces of: bed railings; blood glucose meters; cabinets; carts; chairs; counters; dental unit instrument trays; exam tables; gurneys; isolettes; IV poles; stethoscopes; stretchers; tables; telephones; toilet seats; and hard non-porous outside surfaces of: amalgamators and dental curing lights; diagnostic equipment; patient monitoring equipment; patient support and delivery equipment.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**DIRECTIONS FOR USE** It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**To Dispense Wipes:** Remove lid and discard inner seal from canister. Find center of wipe roll, remove first wipe for use, twist corner of next wipe into a point and thread through the hole in the canister lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out at an angle. When not in use keep center cap of lid closed to prevent moisture loss.

**TO DISINFECT AND DEODORIZE:** To disinfect nonfood contact surfaces only: Unfold a clean wipe and thoroughly wet surface. Allow treated surface to remain wet for a full two (2) minutes. Let air dry. For heavily soiled surfaces, use a wipe to pre-clean prior to disinfecting. These directions also apply to *Mycobacterium bovis* BCG (Tuberculosis) at 68°F (20°C).



(01)10310819000571

## \*SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS

**Personal protection:** When using this product, wear disposable protective gloves, protective gowns, masks, and eye coverings when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids. **Cleaning procedure:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold and use first germicidal wipe to remove heavy soil. **Disposal of infectious materials:** Used wipe, blood and other body fluids should be disposed of according to local regulation for infectious waste disposal. **Contact time:** Use second germicidal wipe to thoroughly wet surface. Allow to remain wet two (2) minutes, let air dry.

### PRECAUTIONARY STATEMENTS Hazards to Humans and Domestic Animals

**WARNING.** Causes substantial but temporary eye damage. Do not get in eyes or on clothing. Avoid contact with skin. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using restroom. Remove and wash contaminated clothing before reuse.

**FIRST AID** Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes.

**PHYSICAL OR CHEMICAL HAZARD** Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

Test wipe on small inconspicuous area first.

**STORAGE AND DISPOSAL** Do not contaminate water, food, or feed by storage and disposal. Storage: Do not store near heat or open flame. When not in use keep center cap of lid closed to prevent moisture loss. **Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Container Disposal:** Nonrefillable container. Do not reuse or refill this container. Offer for recycling. If recycling is not available, put in trash collection.

**EPA REG. NO. 9480-4**

**EPA EST. NO. A= 9480-NY-1, C= 72956-AR-1, D= 8251-WI-4**

**Alpha character will precede batch code on product**

Manufactured by:  
Professional Disposables International, Inc.  
Two Nice-Pak Park  
Orangeburg, NY 10962-1376  
1-800-999-6423  
Made in USA with domestic and imported materials

**NOT FOR USE ON SKIN FOR USE ON HARD, NON POROUS SURFACES ONLY**

# Super

## LARGE WIPE

### WARNING KEEP OUT OF REACH OF CHILDREN

See outer carton for additional  
precautionary statement.

**Directions for use** It is a violation of Federal law to use this product in a manner inconsistent with its labeling. Tear open packet. Unfold wipe and use as directed. Discard after single use.

**TO DISINFECT AND DEODORIZE:** To disinfect nonfood contact surfaces only: Unfold a clean wipe and thoroughly wet surface. Allow treated surface to remain wet for a full two (2) minutes. For heavily soiled surfaces, use a wipe to pre-clean prior to disinfecting. Let air dry. These directions also apply to *Mycobacterium bovis* BCG (Tuberculosis) at 68°F (20°C).

**CONTAINS: 1 LARGE WIPE** 5 x 8 in (12.7 x 20.3 cm) Net Wt. 0.18 oz. (5.1 g)

ACCEPTED  
FOR REGISTRATION

**June 8, 2015**

New York State Department  
of Environmental Conservation  
Division of Materials Management  
Pesticide Product Registration

## Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

### ACTIVE INGREDIENTS:

n-Alkyl (68% C<sub>12</sub>, 32% C<sub>14</sub>) dimethyl  
ethylbenzyl ammonium chlorides ..... 0.25%  
n-Alkyl (60% C<sub>14</sub>, 30% C<sub>16</sub>, 5% C<sub>12</sub>, 5% C<sub>18</sub>)  
dimethyl benzyl ammonium chlorides..... 0.25%  
Isopropyl Alcohol ..... 55.00%

### OTHER INGREDIENTS

TOTAL..... 100.00%

*Does not include the weight of the wipe*

**+PDI**

**PRECAUTIONARY STATEMENTS** Hazards to Humans and Domestic Animals  
**WARNING.** Causes substantial but temporary eye damage. Do not get in eyes or on clothing. Avoid contact with skin. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using restroom. Remove and wash contaminated clothing before reuse.

**FIRST AID** Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. **If in eyes:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes.

**NOT A SKIN WIPE  
FOR USE ON HARD NON POROUS SURFACES ONLY**

**EPA REG. NO.: 9480-4 EPA EST. NO.: 9480-NY-1**

Manufactured by:  
Professional Disposables International, Inc.  
Orangeburg, NY 10962-1376 • USA  
1-800-999-6423



2JH040

<sup>†</sup>SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) OF SURFACES/ OBJECTS SOILED WITH BLOOD/BODY FLUIDS

**PERSONAL PROTECTION:** When using this product, wear disposable protective gloves, protective gowns, masks, and eye coverings when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids.

**CLEANING PROCEDURE:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold and use first germicidal wipe to remove heavy soil.

**DISPOSAL OF INFECTIOUS MATERIALS:** Used wipe, blood and other body fluids should be disposed of according to local regulation for infectious waste disposal.

**CONTACT TIME:** Use second germicidal wipe to thoroughly wet surface. Allow to remain wet two (2) minutes, let air dry.

**PRECAUTIONARY STATEMENTS**  
Hazards to Humans and Domestic Animals

**WARNING.** Causes substantial but temporary eye damage. Do not get in eyes or on clothing. Avoid contact with skin. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using restroom. Remove and wash contaminated clothing before reuse.

**FIRST AID** Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes.

**PHYSICAL OR CHEMICAL HAZARD** Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

Test wipe on small inconspicuous area first.

**STORAGE AND DISPOSAL**  
Do not contaminate water, food, or feed by storage and disposal. **Storage:** Do not store near heat or open flame. **Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Packet Disposal:** Put empty packet in trash collection.



## Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

### ACTIVE INGREDIENTS:

n-Alkyl (68% C<sub>12</sub>, 32% C<sub>14</sub>) dimethyl ethylbenzyl ammonium chlorides.....0.25%

n-Alkyl (60% C<sub>14</sub>, 30% C<sub>16</sub>, 5% C<sub>12</sub>, 5% C<sub>18</sub>) dimethyl benzyl ammonium chlorides.....0.25%

Isopropyl Alcohol.....55.00%

**OTHER INGREDIENTS** .....44.50%

**TOTAL** .....100.00%

*Does not include the weight of the wipe.*

## WARNING KEEP OUT OF REACH OF CHILDREN

See back panel for additional precautionary statements

**CONTAINS: 50 Wipes**  
5.0 x 8.0 in (12.7 x 20.3 cm)  
Net Wt. 9.0 oz (257g)



## ACCEPTED FOR REGISTRATION

June 8, 2015

New York State Department  
of Environmental Conservation  
Division of Materials Management  
Pesticide Product Registration

5JH04095

REORDER NO. H04082

### Effective against:

**BACTERIA** • *Bordetella pertussis*  
• *Burkholderia cepacia* • *Campylobacter jejuni*  
• *Escherichia coli* [E. coli] • *Escherichia coli* O157:H7 • *Klebsiella pneumoniae*  
• *Pseudomonas aeruginosa* • *Salmonella enterica* • *Staphylococcus aureus*  
**MULTI-DRUG RESISTANT BACTERIA**  
• *Acinetobacter baumannii* • *Enterobacter cloacae* NDM-1 Positive • ESBL Producing *Escherichia coli* [E. coli] • *Klebsiella pneumoniae* - KPC-2 positive, ST258  
• Methicillin Resistant *Staphylococcus aureus* (MRSA) • Vancomycin Resistant *Enterococcus faecalis* (VRE)

**VIRUSES\*** • \*Adenovirus Type 5 • \*Herpes Simplex Virus Type 2 • \*Human Coronavirus Strain 229E • \*Influenza A (H1N1) Virus  
• \*Influenza A/Hong Kong • \*Respiratory Syncytial Virus (RSV) • \*Rhinovirus  
• \*Rotavirus Strain WA • \*Vaccinia Virus  
• Kills Pandemic 2009 H1N1 Influenza A virus  
**BLOODBORNE PATHOGENS** • \*Hepatitis B Virus (HBV) - Duck HBV† • \*Hepatitis C Virus (Human)(HCV) - Bovine Diarrhea Virus†  
• \*HIV-1(AIDS virus)†

**TB** • *Mycobacterium bovis* - BCG (Tuberculosis)  
**PATHOGENIC FUNGI** • *Candida albicans*

**NOT FOR USE ON SKIN**  
**FOR USE ON HARD NON POROUS SURFACES ONLY**

**EPA REG. NO.: 9480-4**  
**EPA EST. NO.: 9480-NY-1**

Manufactured by:  
Professional Disposables International, Inc.  
Two Nice-Pak Park, Orangeburg, NY 10962-1376 • USA  
1-800-999-6423

## Super

LARGE WIPE

## Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

### BACTERICIDAL, TUBERCULOCIDAL, VIRUCIDAL\* IN 2 MINUTES SUITABLE FOR USE ON EQUIPMENT REQUIRING ALCOHOL-BASED PRODUCTS

<sup>†</sup>KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature (68°-77°F) in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

**AREAS OF USE:**  
**Hospital and Healthcare Settings:** Ambulatory Surgical Centers (ASC), Clinics, Dental Offices, Dialysis Clinics, Home Health Care, Hospices, Hospitals, Laboratory, Nursing homes, Physical therapy, Physician's offices, Radiology, Rehabilitation, Transport vehicles.

**CRITICAL CARE AREAS:** CCU, Emergency Rooms, ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU), Surgery and Surgical Intensive Care Units (SICU).

**Hospital, Healthcare, and Critical Care Use Sites:** May be used on hard non-porous surfaces of: bed railings, blood glucose meters, cabinets, carts, chairs, counters, dental unit instrument trays, exam tables, gurneys, isolettes, IV poles, stethoscopes, stretchers, tables, telephones,

toilet seats; and hard non-porous outside surfaces of: amalgamators and dental curing lights, diagnostic equipment, patient monitoring equipment, patient support and delivery equipment.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**DIRECTIONS FOR USE** It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**BOX:** Open along perforations as indicated. Remove packets as needed. **PACKET:** Tear open packet. Unfold and use as directed. Discard after single use.

**TO DISINFECT AND DEODORIZE:** To disinfect nonfood contact surfaces only: Unfold a clean wipe and thoroughly wet surface. Allow treated surface to remain wet for a full two (2) minutes. Let air dry. For heavily soiled surfaces, use a wipe to pre-clean prior to disinfecting. These directions also apply to *Mycobacterium bovis* BCG (Tuberculosis) at 68°F (20°C).

# Super

## EXTRA LARGE WIPE

# Sani-Cloth®

### GERMICIDAL DISPOSABLE WIPE

#### ACTIVE INGREDIENTS:

n-Alkyl [68% C <sub>12</sub> , 32% C <sub>14</sub> ] dimethyl ethylbenzyl ammonium chlorides.....	0.25%
n-Alkyl [60% C <sub>12</sub> , 30% C <sub>14</sub> , 5% C <sub>12</sub> , 5% C <sub>14</sub> ] dimethyl benzyl ammonium chlorides.....	0.25%
Isopropyl Alcohol.....	55.00%
<b>OTHER INGREDIENTS</b> .....	44.50%
<b>TOTAL</b> .....	100.00%

*Does not include the weight of the wipe*

**KEEP OUT OF REACH OF CHILDREN**

**WARNING** See back panel for additional precautionary statements

\*KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature [68°-77°F] in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) [associated with AIDS], Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

**Areas of Use: Hospital and Healthcare Settings:** Ambulatory Surgical Centers (ASC), Clinics, Dental Offices, Dialysis Clinics, Home Health Care, Hospices, Hospitals, Laboratory, Nursing homes, Physical therapy, Physician's offices, Radiology, Rehabilitation, Transport vehicles.

**Critical Care Areas:** CCU, Emergency Rooms, ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU), Surgery and Surgical Intensive Care Units (SICU). **Hospital, Healthcare, and Critical Care Use**

**Sites:** May be used on hard non-porous surfaces of: bed railings, blood glucose meters, cabinets, carts, chairs, counters, dental unit instrument trays, exam tables, gurneys, isolettes, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats; and hard non-porous outside surfaces of: amalgamators and dental curing lights, diagnostic equipment, patient monitoring equipment, patient support and delivery equipment.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**CONTAINS: 1 EXTRA LARGE WIPE 11.5 x 11.75 in (29.2 x 29.8 cm) Net Wt. 0.5 oz (16g)**

REORDER NO. U872  
ACCEPTED  
FOR REGISTRATION  
**June 8, 2015**

New York State Department  
of Environmental Conservation  
Division of Materials Management  
Pesticide Product Registration

**Packet:** Tear open packet. Unfold wipe and use as directed. Discard after single use.

**Directions for use:** It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**TO DISINFECT AND DEODORIZE:** To disinfect nonfood contact surfaces only: Unfold a clean wipe and thoroughly wet surface. Allow treated surface to remain wet for a full two (2) minutes. Let air dry. For heavily soiled surfaces, use a wipe to pre-clean prior to disinfecting. These directions also apply to *Mycobacterium bovis* BCG (Tuberculosis) at 68°F (20°C).

**\*SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS**

**Personal protection:** When using this product, wear disposable protective gloves, protective gowns, masks, and eye coverings when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids.

**Cleaning procedure:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold and use first germicidal wipe to remove heavy soil.

**Disposal of infectious materials:** Used wipe, blood and other body fluids should be disposed of according to local regulation for infectious waste disposal.

**Contact time:** Use second germicidal wipe to thoroughly wet surface. Allow to remain wet two (2) minutes, let air dry.

**PRECAUTIONARY STATEMENTS** Hazards to Humans and Domestic Animals

**WARNING:** Causes substantial but temporary eye damage. Do not get in eyes or on clothing. Avoid contact with skin. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using restroom. Remove and wash contaminated clothing before reuse.

**FIRST AID** Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes.

**PHYSICAL OR CHEMICAL HAZARD** Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

Test wipe on small inconspicuous area first.

**STORAGE AND DISPOSAL:** Do not contaminate water, food, or feed by storage and disposal. **Storage:** Do not store near heat or open flame. **Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Container Disposal:** Put empty packet in trash collection.

**NOT FOR USE ON SKIN**

**FOR USE ON HARD,  
NON POROUS SURFACES ONLY**

**EPA REG. NO.: 9480-4**

**EPA EST. NO.: 9480-NY-1**

Manufactured by:  
Professional Disposables  
International, Inc.

Orangeburg, NY 10962-1376 • USA  
1-800-999-6423



(01)10310819000601

2NU872

**+PDI**

2NU872\_2015\_05\_12.pdf





GERMICIDAL DISPOSABLE WIPE

Sani-Cloth®

EXTRA LARGE  
WIPE

Super

5MU87250

Super

Sani-Cloth®  
GERMICIDAL DISPOSABLE WIPE

EXTRA LARGE  
WIPE

BACTERICIDAL, TUBERCULOCIDAL, VIRUCIDAL\* IN 2 MINUTES

SUITABLE FOR USE ON EQUIPMENT REQUIRING  
ALCOHOL-BASED PRODUCTS

\*KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON THE PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS AND PROBES PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in two (2) minutes at room temperature (68°-77°F) in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

#### AREAS OF USE:

**Hospital and Healthcare Settings:** Ambulatory Surgical Centers (ASC), Clinics, Dental Offices, Dialysis Clinics, Home Health Care, Hospices, Hospitals, Laboratory, Nursing homes, Physical therapy, Physician's offices, Radiology, Rehabilitation, Transport vehicles.

**CRITICAL CARE AREAS:** CCU, Emergency Rooms, ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU), Surgery and Surgical Intensive Care Units (SICU).

**Hospital, Healthcare, and Critical Care Use Sites:** May be used on hard non-porous surfaces of: bed railings, blood glucose meters, cabinets, carts, chairs, counters, dental unit instrument trays, exam tables, gurneys, isolettes, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats; and hard non-porous outside surfaces of: amalgamators and dental curing lights, diagnostic equipment, patient monitoring equipment, patient support and delivery equipment. This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**Box:** Open along perforations as indicated. Remove packets as needed.

**DIRECTIONS FOR USE:** It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**TO DISINFECT AND DEODORIZE:** To disinfect nonfood contact surfaces only: Unfold a clean wipe and thoroughly wet surface. Allow treated surface to remain wet for a full two (2) minutes. Let air dry. For heavily soiled surfaces, use a wipe to pre-clean prior to disinfecting. These directions also apply to *Mycobacterium bovis* BCG (Tuberculosis) at 68°F (20°C).

\*SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS

**Personal protection:** When using this product, wear disposable protective gloves, protective gowns, masks, and eye coverings when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids.

**Cleaning procedure:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal wipe. Open, unfold and use first germicidal wipe to remove heavy soil.

**Disposal of infectious materials:** Used wipe, blood and other body fluids should be disposed of according to local regulation for infectious waste disposal.

**Contact time:** Use second germicidal wipe to thoroughly wet surface. Allow to remain wet two (2) minutes, let air dry.

**PRECAUTIONARY STATEMENTS** Hazards to Humans and Domestic Animals

**WARNING.** Causes substantial but temporary eye damage. Do not get in eyes or on clothing. Avoid contact with skin. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using restroom. Remove and wash contaminated clothing before reuse.

**FIRST AID** Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. If in eyes: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15-20 minutes.

**PHYSICAL OR CHEMICAL HAZARD** Combustible. Do not use or store near heat or open flame. Do not use on natural marble, windows, unpainted wood, brass, clear plastic or colored grout.

Test wipe on small inconspicuous area first.

#### STORAGE AND DISPOSAL

Do not contaminate water, food, or feed by storage and disposal.

**Storage:** Do not store near heat or open flame.

**Towelette Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet.

**Packet Disposal:** Put empty packet in trash collection.

ACCEPTED  
FOR REGISTRATION

June 8, 2015

New York State Department  
of Environmental Conservation  
Division of Materials Management  
Pesticide Product Registration

Super

Sani-Cloth®  
GERMICIDAL DISPOSABLE WIPE

EXTRA LARGE  
WIPE



BACTERICIDAL, TUBERCULOCIDAL, VIRUCIDAL\* IN 2 MINUTES

SUITABLE FOR USE ON EQUIPMENT REQUIRING  
ALCOHOL-BASED PRODUCTS

#### ACTIVE INGREDIENTS:

n-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> ) dimethyl ethylbenzyl ammonium chlorides .....	0.25%
n-Alkyl (60% C <sub>14</sub> , 30% C <sub>16</sub> , 5% C <sub>12</sub> , 5% C <sub>18</sub> ) dimethyl benzyl ammonium chlorides.....	0.25%
Isopropyl Alcohol .....	55.00%
<b>OTHER INGREDIENTS</b> .....	44.50%
<b>TOTAL</b> .....	100.00%

*Does not include the weight of the wipe.*

**KEEP OUT OF REACH OF CHILDREN**  
**WARNING** See back panel for additional  
precautionary statements

TEAR AT PERFORATION

**CONTAINS: 50 Wipes**  
11.5 x 11.75 in (29.2 x 29.8 cm)  
Net Wt. 1.9 lbs (872g)



REORDER NO. U87295

#### Effective against:

**BACTERIA** • *Bordetella pertussis* • *Burkholderia cepacia* • *Campylobacter jejuni* • *Escherichia coli* (*E. coli*) • *Escherichia coli* O157:H7 • *Klebsiella pneumoniae* • *Pseudomonas aeruginosa* • *Salmonella enterica* • *Staphylococcus aureus*

#### MULTI-DRUG RESISTANT BACTERIA

• *Acinetobacter baumannii* • *Enterobacter cloacae* NDM-1 Positive • ESBL Producing *Escherichia coli* (*E. coli*) • *Klebsiella pneumoniae* - KPC-2 positive, ST258 (multi-drug and carbapenem resistant) • Methicillin Resistant *Staphylococcus aureus* (MRSA) • Vancomycin Resistant *Enterococcus faecalis* (VRE)

**VIRUSES\*** • \*Adenovirus Type 5 • \*Herpes Simplex Virus Type 2 • \*Human Coronavirus Strain 229E • \*Influenza A (H1N1) Virus • \*Influenza A/Hong Kong • \*Respiratory Syncytial Virus (RSV) • \*Rhinovirus • \*Rotavirus Strain WA • \*Vaccinia Virus • Kills Pandemic 2009 H1N1 Influenza A virus

#### BLOODBORNE PATHOGENS

• \*Hepatitis B Virus (HBV) - Duck HBV\* • \*Hepatitis C Virus (Human)(HCV) - Bovine Diarrhea Virus\* • \*HIV-1(AIDS virus)\*

**TB** • *Mycobacterium bovis* - BCG (Tuberculosis)

**PATHOGENIC FUNGI** • *Candida albicans*

NOT FOR USE ON SKIN  
FOR USE ON HARD NON POROUS SURFACES ONLY

EPA REG. NO.: 9480-4  
EPA EST. NO.: 9480-NY-1

Manufactured By:  
Professional Disposables International, Inc.  
Two Nice-Pak Park, Orangeburg, NY 10962-1376 • USA  
1-800-999-6423

